
Pulse Oximeter User Manual

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Dear customer,

Thank you for purchasing our product. Please read the following instructions carefully before first use and keep this user manual for future reference. Pay particular attention to the safety instructions. If you have any questions or comments about the device, please contact the customer line.

www.alza.co.uk/kontakt

+44 (0)203 514 4411

Importer

Alza.cz a.s., Jankovcova 1522/53, Holešovice, 170 00 Praha 7, www.alza.cz

Safety

Warnings

DO NOT use the oximeter while the patient is being scanned by MRI.

It is recommended that the oximeter should be inspected before use, when there is obvious damage, stop using the device.

Discomfort or pain may appear if using the device ceaselessly. Especially for microcirculation barrier patients, it recommended that the oximeter should not be used on the same finger more than 30 minutes.

For some patients, need a more careful examination of the measurement site, the oximeter shall not be placed in edema or fragile organization E. The device is not used for neonate and infant.

The device is just one of clinical auxiliary equipment, other means of clinical diagnosis and clinical response factors determine the patient's condition must be consolidated under the guidance of a doctor.

It is not recommended to use the oximeter in high frequency environment such as electrosurgical equipment.

Prevent children from swallowing the oximeter or its accessories children must be, accompanied by adult guardian using products.

Please follow local ordinances and recycling instructions regarding disposal or recycling of the device and batteries.

Cautions

The equipment is designed to measure the percentage of arterial oxygen saturation of functional hemoglobin. Factors that may degrade pulse oximeter performance or affect the accuracy of the measurement include the following:

Do not apply the pulse oximeter on the same arm as a blood pressure cuff, arterial catheter, or infusion line(s)(IVs).

- Excessive light, such as sunlight or direct home lighting.
- Moisture in the equipment.
- Finger is outside recommended size range.
- Poor pulse quality.
- Venous pulsations.
- Anemia or low hemoglobin concentrations.
- Cardio green and other intravascular dyes.
- Carboxyhemoglobin.
- Methemoglobin.
- Dysfunctional hemoglobin.

- Artificial nails or fingernail polish.

Please read the measure value when the numbers or waveform on screen is equably and steady going, the measure value is optimal value.

The light (the infrared is invisible) emitted from the equipment is harmful to the eyes, so the user and the maintenance man should not look at the light source.

The principle

The principle of pulse oximeter is based on the red and infrared light absorption characteristics of oxygenated and deoxygenated hemoglobin. Oxygenated hemoglobin absorbs more infrared light and allows more red light to pass through (Figure 1). Deoxygenated (or reduced) hemoglobin absorbs more red light and allows more infrared light to pass through. Red light is in the 600-750 nm wavelength light band. Infrared light is in the 850-1000 nm wavelength light band.

Pulse oximeter sensors have red and infrared low voltage light emitting diodes (LEDs) which serve as light sources. The emitted light is transmitted through the tissue, then detected by the photodetector and sent to the microprocessor of the pulse oximeter. All constituents of the human body, venous and arterial blood, and tissue absorb light (Figure 2). The pulsating of arterial blood results in changes in the absorption to added hemoglobin (Hb) and oxygenated hemoglobin (HbO₂) in the path of the light. Since HbO₂ and Hb absorb light to varying degrees, this varying absorption is translated into plethysmographic waveforms at both red and infrared wavelengths. The relationship of red and infrared plethysmographic signal amplitude can be directly related to arterial oxygen saturation.


Figure 1 Absorption spectra of HbO₂ and Hb


Figure 2 Oximeter's schematic illustration for red and infrared wavelengths

Intended Use

The pulse oximeter is intended for measure oxygen saturation and pulse rate of adults or children in families, hospitals, oxygen bars, social medical care institutions, you can use this equipment for measurement before or after sports.

Product Introduction and Operation Guide

Front View


Figure 3 Front View of YM101/YM102/YM103

Operation Method

Open the battery cover, and put the two AAA batteries into the battery compartment incorrect polarities, then replace the cover;

Press the bottom of the equipment and open the probe, then insert one finger into the probe;

Press the button to turn the equipment on, and the measure interface will appear;

After about 8 seconds, the measurement result can be read directly from the display screen;

Before reading the parameters, make sure that stable numbers of the pulse oximeter interface has sustained more than 4 second;

The equipment will turned off automatically within 8 seconds when the finger left the probe.

Battery Installation

- Put the two AAA batteries into battery compartment in correct polarities (Figure4).
- Push the battery cover horizontally along the arrow shown as right.

WARNINGS:

- Battery polarities should be correctly installed, otherwise, damage may be caused to the equipment.
- Please remove the batteries if the equipment will not use for a long time.

Lanyard installation

Pass the thinner end of the lanyard through the hanging hole;

Pass the thicker end of the lanyard through the thinner end and tighten the lanyard (Figure5).


Figure 4 Battery Installation


Figure 5 Lanyard Installation

Attention for Operation

Before use check and confirm that the people or finger size were applicable;

Before use check and confirm that the environment should be non-combustible material, as well as to avoid high or low temperature and humidity, but also need to pay attention to the following:

To avoid glare and direct sunlight exposure;

To avoid radiation infrared or ultraviolet radiation;

Avoid contact with the organic solvent, mist, dust, corrosive gases;

The equipment should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection;

The equipment may not work normally on microcirculation barrier patients, Warm or rub the finger, or re-position the equipment could improve the measurement. The ray between photo detector and light emitting diode should across patient's arteriole. F. The patient should not use enamel or other makeup;

Avoid inserting a wet finger into the probe.

Notes:

- The user should fully insert the finger into the probe,

- It is recommended to let the LED light shine directly on the nail (Figure 6);
- Don't shake the finger and try to keep the patient still during the measurement.


Figure 6 Finger Placement Diagram

Functions and Menu Operation

The button operation rules

Long-press functions include entering menus, activating item's submenu, confirming setting values, and exiting item's submenu; short-press functions are polling menu items and viewing the setting values of items. It should be noted that long-press means pressing the key for about 2 seconds, and short-press means pressing the key for less than 0.5 second.

Switch display direction

Short press to switch the display direction.




Menu Operation Active the menu

After the oximeter is turned on, long-press the power button to activate the menu, then short press the button to view the setting values of each item. If the user wants to change the setting value of the item, long-press to enter the item's submenu, the parameter value starts to flash, short-press to traverse the parameter value until the parameter value required by the user is selected, long press to confirm and exit the submenu.



Item1. Setup the LED display brightness

The first item is to setup the display brightness. Long-press the button to select a brightness level ranging from 1 to 3. The greater the value, the greater brightness of the display.


Item2.Setup the SpO2 Alarm Limits

The second item is to setup the SpO2 alarm limits. For example: When Spo2 High limit is set to 96, an alarm will be issued when the spo2 value is higher than 92, and when Spo2 low limit is set to 94, an alarm will be issued when the spo2 value is lower than 94.


Item3.Setup the PR Alarm Limits

The third item is to setup the PR alarm limits. For example: When PR High limit is set to 130, an alarm will be issued when the PR value is higher than 130, and when PR low limit is set to 50, an alarm will be issued when the PR value is lower than 50.

Item4.Turn Sound On/Off

The forth item is short-press to turn sound on /off.

Specifications

Classification

Type of protection against electric shock	II(Internally powered equipment)
Degree of protection against electric shock	Type BF-Applied part (non-defibrillation proof)
Operating mode	Spot checking
Degree of protection against hazards of explosion	IP22

Power Requirements

Specification of alkaline batteries	Two AAA (LR03)
Operating current	25-50mA

Physical Specifications

Width*Height*Depth	57×30×31 mm
Weight	28g (Bare machine)

Measurement Specifications

Spo2 declared accuracy

Range	70%~100%: ±2digits 0% ~ 69%: unspecified
Resolution	1%
SpO2 Display Range	30%~99%
S 2A PO accuracy under conditions of low perfusion	±3% when PI is 0.075%

PR declared accuracy

Range	25~250: ±3digits
Resolution	1 bpm
PR Accuracy under conditions of low perfusion	±3bpm when Pi is 0.075%

Environmental Specifications

Temperature

Operating	+50° to +104°F / +10° to +40°C
Storage/Transportation	-4° to +140°F / -20° to +60°C

Humidity

Operating	15~95%, noncondensing
Storage/Transportation	10~95%, noncondensing

Atmosphere Pressure

Operating	70~106kpa
Storage/Transportation	50~107.4kpa

Display

Model	Display type	Display content
YM101/YM102/YM103	1.5" LED Display; YM101: Red YM102: Green YM103: White	SpO2%, Pulse Rate, PI%, Bar Graph Battery Indicator, Heartbeat & Buzzer Icon

Notes:

The claim for oxygen saturation accuracy should be supported by clinical studies covering the entire claimed range, the fraction of inspired oxygen (FiO2) delivered to test subjects is varied to achieve a series of targeted steady-state saturation periods over the specified SpO2 accuracy range (e.g. 70 % to 100 %), then the SpO2 accuracy is calculated by comparing SpO2 readings of the pulse oximeter to the values of SaO2 determined with a Co-Oximeter.

The clinical trial included 11 subjects, including 6 males and 5 females, with an age range of 18 to 46 years, the subjects skin color included dark black, medium black, light color and white.

The accuracy of SPO2 and PR of the oximeter under weak perfusion has been verified using the signal of Fluke Index2 SpO2 SIMULATOR. The SPO2 and PR values of the oximeter differ from those under weak signal conditions and are compared to the input signals of known SPO2 and PR.

Maintenance, Cleaning

Maintenance

The equipment's design life expectancy is about 2 years, keep your equipment and accessories free of dust and dirt, and follow these rules: A. Please clean the equipment before use according to chapter 6.2; Remove the batteries inside the battery cassette if the equipment will not be operated for a long time; B. Replace the batteries in time when the battery voltage indicate lamps were empty; C. It is recommended that the equipment should be kept in a dry environment with no corrosive gases and good ventilation anytime. The moisture and high-light environments will affect its lifetime and even might damage the equipment. D. It is best to preserve the product in a place where the temperature is between -20 to 60°C and the relative humidity is less than 95%. E. The packed equipment can be transported by ordinary conveyance. The equipment is not transported mixed with toxic, harmful, corrosive materials.

WARNING

No modification of this equipment is allowed.

Cleaning

Your equipment should be cleaned on a regular basis. If there is heavy pollution or lots of dust and sand in your place, the equipment should be cleaned more frequently. Before cleaning the sand in your place, the equipment should be cleaned more frequently. Before cleaning the cleaning, agents are:

- a) Mild soap (diluted).
- b) Ethanol (70%).

To clean your equipment, follow these rules:

- a Shut down the pulse oximeter;
- b Clean the display screen using a soft, clean cloth dampened with a glass cleaner;
- c Clean the exterior surface of the equipment and probe using a soft cloth dampened with the cleaner;
- d Wipe off all the cleaning solution with a dry cloth after cleaning if necessary;
- e Dry your equipment in a ventilated, cool place.

To avoid damage to the equipment, follow these rules:

CAUTIONS

Always dilute according to the manufacturer's instructions or use lowest possible concentration. Do not immerse part of the equipment in the liquid. Do not pour liquid onto the equipment or accessories. Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners). If you spill liquid onto the equipment, contact us or your service personnel.

Disposal

Dispose of the pulse oximeter in accordance with local environment and waste disposal laws and regulations.

Accessories

One lanyard.

Two AAA batteries (Optional).



One user manual.

One certificate card.

Troubleshooting

Trouble	Possible Reason	Solutions
The device cannot be turned on	The battery is drained away or almost drained away.	Please replace battery
	The battery installation is incorrect.	Install the battery over again
	The device works abnormally.	Please contact the product distributor
The display is off suddenly	The device was set to shut down automatically in 8 seconds when there is no correct physiological signals	Normal
	The battery is almost drained away	Please replace the batteries
The Spo2 and PR are not displayed stably.	The finger is not inserted deep enough	Replace the finger and try again
	The finger is shaking, or the body is moving	Try to keep still
	Not used in the work environment required by this manual	Please use in normal working environment
	The device works abnormally	Please contact the product distributor
The Spo2 and PR are not displayed normally	The finger size is too big or small	Select the suitable size finger to measure
	.Excessive ambient light	Avoid the excessive ambient light irradiation
	User's blood perfusion is very low	Warm the finger and try again

Symbols and Definition

Symbol	Definition	Symbol	Definition
%SpO2	The Pulse Oxygen Saturation (%)		The device has no Alarm System
PRbpm	The Pulse Rate (beat per Minute)		Date of Manufacturer
PI%	The blood perfusion index (%)		Manufacturer
	This item is compliant with Medical Device Directive 93/42/EEC	SN	Serial Number
IP22	Degrees of protection provided by enclosures per IEC60529		Keep away from rain
	Battery positive electrode		Battery cathode electrode
	BF type applied part		This way up
	Caution, consult accompanying documents		Stacking limit by number
	Power button		Humidity limitation
	Low voltage indicator		Temperature limitation
	Atmospheric pressure limitation		

Warranty Conditions

A new product purchased in the Alza.cz sales network is guaranteed for 2 years. If you need repair or other services during the warranty period, contact the product seller directly, you must provide the original proof of purchase with the date of purchase.

The following are considered to be a conflict with the warranty conditions, for which the claimed claim may not be recognized:

- Using the product for any purpose other than that for which the product is intended or failing to follow the instructions for maintenance, operation and service of the product.
- Damage to the product by a natural disaster, the intervention of an unauthorized person or mechanically through the fault of the buyer (e.g. during transport, cleaning by inappropriate means, etc.).
- Natural wear and aging of consumables or components during use (such as batteries, etc.).
- Exposure to adverse external influences, such as sunlight and other radiation or electromagnetic fields, fluid intrusion, object intrusion, mains overvoltage, electrostatic discharge voltage (including lightning), faulty supply or input voltage and inappropriate polarity of this voltage, chemical processes such as used power supplies, etc.
- If anyone has made modifications, modifications, alterations to the design or adaptation to change or extend the functions of the product compared to the purchased design or use of non-original components.

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