English

90mm

65mm

User's manual for electric toothbrush

As recognized by a large number of customers, we commit to achieving scientific and universal oral health for the communities that we base. Please read the product-specific details thoroughly for better services.

*** IMPORTANT ***

Before operating this toothbrush, please read these instructions carefully. Any comments or suggestions are welcome!

Brush head compartment

Brush head

PTO shaft

On/Off

Brush mode button

Brushing mode display

Charge level display

USB-DC Cable

Name

Electric Toothbrush

Model

J20C

Waterproof level

IPX7

Rated voltage and power

DC3.7V-0.5W

Charging time

4 Hours

Standards

GB4706.1-2005, GB4706.59-2008

Accessories

Handle

*1

Brush head

*2

User manual

*1

USB Cable

*1

Important information on charging!

1) Verify voltage and outlets are safe and normal. Please use correct power adapter when charging.

2) Make sure that charging cable is dry, charging port free from water stains and damage, etc. Connect charger to a power outlet when it is safe.

3) If this is your first time using the product, you should test it. Fully charge the battery when it is low or empty.

Toothbrush Battery Indicator

Toothbrush Status

Light Indicator

Brushing

Light off

Low Battery

Flashing Red

Charging

Solid Red

Full Charged

From Red to Solid Green

Notes:

1) Please be attention with the battery light status when using the toothbrush.

2) Charge it timely in case the battery is empty or runs low.

If the toothbrush is not used for a long time, please pre-charge before operating it.Safe and correct charging can increase the lifespan of the product.

Warranty and support

If you need information or support, please visit www.alza.co.uk/contact or call us at +44 (0)203 514 4411. **Importer:** Alza.cz a.s., Jankovcova 1522/53, Holešovice, 170 00 Praha 7, www.alza.cz.

A new product purchased in the Alza.cz sales network is guaranteed for 2 years. If you need repair or other services during the warranty period, contact the product seller directly, you must provide the original proof of purchase with the date of purchase.

The following are considered to be a conflict with the warranty conditions, for which the claimed claim may not be recognized:

1) Using the product for any purpose other than that for which the product is intended or failing to follow the instructions for maintenance, operation, and service of the product.

2) Damage to the product by a natural disaster, the intervention of an unauthorized person or mechanically through the fault of the buyer (e.g. during transport, cleaning by inappropriate means, etc.).

3) Natural wear and aging of consumables or components during use (such as batteries, etc.).

4) Exposure to adverse external influences, such as sunlight and other radiation or electromagnetic fields, fluid intrusion, object intrusion, mains overvoltage, electrostatic discharge voltage (including lightning), faulty supply or input voltage and inappropriate polarity of this voltage, chemical processes such as used power supplies, etc.

5) If anyone has made modifications, modifications, alterations to the design or adaptation to charge or extend the functions of the product compared to the purchased design or use of non-original components.

This equipment is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU, 2014/30/EU and Directive 2011/65/EU amended (EU) 2015/863.

This product must not be disposed of as normal household waste in accordance with the EU Directive on Waste Electrical and Electronic Equipment (WEEE - 2012/19 / EU).

CE

Czech

90mm

65mm

User's manual for electric toothbrush

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