BerlingerHaus®
PROFESSIONAL PRECISION

BH-9100, BH-9101, BH-9102, BH-9108, BH-9220,
BH-9223, BH-9356, BH-9358, BH-9527, BH-9528, BH-9684

DIGITAL BATHROOM BODYFAT SCALE

INSTRUCTION MANUAL


GB, RO, SK, CZ, FR, PL, ES, NL, LT, DE, IT, HU, RU

USER MANUAL

GENERAL SAFETY INSTRUCTIONS

LCD
display

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing. If you give this device to other people, please also pass on the operating instructions.

- The device is designed exclusively for private and not for commercial use.
- Only use the device for the intended purpose.
- Do not use a damaged device.
- Keep away from heat, direct sunlight, damp (never immerse in liquids) and sharp edges.
- Do not try to repair the appliance on your own. Always contact an authorized technician.
- Pay careful attention to the following "Special Safety Instructions".

Children

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.

WARNING!

- Do not allow small children to play with the foil as there is a danger of suffocation!
- Children should be supervised at all times in order to ensure that they do not play with the device.

Special Safety Instructions For This Machine

WARNING:

If the glass surface is moist or covered with water droplets, please wipe the surface dry as there is a danger of slipping.

CAUTION:

- Handle the balance with care, as it is a precise measuring instrument. Do not drop it or jump on it.
- Do not load the appliance with more than 180 kg to avoid damaging the device.
- If the scales are carried or not going to be used, please store them only in a horizontal position.
- The body fat analysis scale is not suitable for the following persons:
 - Children under 10 years and adults over 99 years
 - Persons with symptoms of fever, edema, or osteoporosis
 - Persons undergoing treatment for dialysis
 - Persons who use a pacemaker
 - Persons who take cardiovascular medication
 - Pregnant women
 - Sports persons who perform over 5 hours of intensive sporting activities each day
 - Persons with a heartbeat lower than 60 beats/minute

Handling batteries

WARNING:

- Do not expose the batteries to high temperatures or direct sunlight. Never throw batteries into the fire. There is a danger of explosion!
- Keep batteries away from children. They are not a toy!
- Do not open the package of the batteries by bigger force.

WARNING:

- Avoid contact with metallic objects (Rings, nails, screws etc.) There is a danger of short-circuiting!
- As a result of a short-circuit batteries may heat up considerably or even catch fire. This may result in burns.
- For your own safety the battery poles should be covered with adhesive strips during transport.
- If a battery leaks do not rub the liquid into your eyes or mucous membranes. If hands come into contact with the liquid they should be washed thoroughly and eyes rinsed out with clear water; if symptoms persist please consult a doctor.

Inserting the batteries

NOTE:

Where required, the supplied batteries have been protected for transport with plastic foil.

1. Remove this foil before inserting the batteries into the appliance.
2. Open the battery compartment on the back.
3. Insert 2 batteries of type AAA/R03. Please ensure that the polarity is correct (see the look of the battery compartment)!
4. Close the battery compartment again.
5. If the device is not used for a prolonged period of time, please take the batteries out in order to prevent leakage of battery acid.
6. If the message "Lo" appears on the display, replace the batteries.

CAUTION:

- Different battery types or new and used batteries must not be used together.
- Batteries are not to be disposed of together with domestic wastes. Dispose of the used batteries at the relevant collection bin.

Location of controls

1. Contact Surfaces
2. Display
3. DOWN-button ▼
4. UP-button ▲
5. SET-button

Display Symbol	Meaning
Kg (lb/st)	Body Weight
USER	User memory space
♂	Male
♀	Female
cm (ft/in)	Body height
Age	Age
☞	Body fat display
☞	Body water display
☞	Muscle mass display
☞	Bone weight display
%	Percentage-body fat/body water/muscle mass/bone weight
KCAL	Daily kilocalorie requirement
Increase	Weight gain in kg(lb/st) compared with previous measurement
Decrease	Weight loss in kg(lb/st) compared with previous measurement

Technical Specifications

Model:756
Power supply:2 Batteries 1.5 V Type "AAA" / "R03"
Net weight of device:1.3 kg
Memory capacity:10 P0-P9
Mass alignment:from 100 to 220 cm
Age:from 18 to 100 years
Minimum weight:from 3.0 kg
Maximum weight:180 kg
Weight measurement unit:100g increments
Body fat measurement unit:0.1 %
Body fat measuring range:4.0 % - 50.0 %
Body water measuring range:35.0 % - 85.0 %
Muscle mass measuring range:10.0 % - 50.0 %

Subject to technical changes without prior notice!

This device has been tested according to all relevant current CE guidelines, such as electromagnetic compatibility and low voltage directives, and has been constructed in accordance with the latest safety regulations.

Information on Measuring Body Fat

WARNING:

Do not use the scales if you use a pacemaker or another medical implant with electronic components. The function of the implant may be impaired by the measurement.

NOTE:

The main principle of this device is based on the measurement of the electric impedance within the human body. In order to determine the actual body data, the device leads a weak, unnoticeable electric signal through your body.

The signal is generally safe! However, please read our special safety instructions!

The method of measurement is referred to as "bioelectric impedance analysis" (BIA), a factor which is based on the relationship between body weight and associated body water percentages as well as on other biological data (age, gender, body weight).

There are heavy persons with a well-balanced percentage between fat and muscle mass, and there are seemingly slim persons with significantly high body fat values. Fat is a factor that puts serious strain on the body, in particular the cardiovascular system. Therefore, controlling both values (fat and weight) is a prerequisite for health and fitness.

Body Fat

Body fat is an important component of the body. From protecting the joints as a vitamin repository to regulating the body temperature, fat plays a vital role in the body. Thus the goal is not to radically reduce the proportion of body fat, but rather to achieve and maintain a healthy balance between muscle and fat.

The Importance of Water for our Health

Approximately fifty-five to sixty percent of the human body consists of water, depending on age and gender in varying amounts.

Water has to fulfill several tasks within the body:

- It is the building block of our cells. That means that all the cells in our body, whether skin cells, gland cells, muscle cells, brain cells or others, are functional only if they contain enough water.
- It is a solvent for important substances in our body.
- It is a means of conveyance for nutrients, the body's own substances and metabolic products.

Reference values for your fitness profile can be seen in the benchmark table. It shows the dependency of body fat data on age and gender. Always pay attention to your body's water balance, which, as well as the fat level, undergoes strong fluctuations throughout the day. Therefore, to obtain consistent data, always measure at the same time of day.

The calculated results of the body fat and water level should only be considered as reference values. They are not recommended to be used for medical purposes. Please always consult a doctor if you have questions about a medical treatment or a diet in order to achieve a normal body weight.

Guideline values for body fat and body water content as a percentage (approximate figures)

Male	10-20	21-42	43-65	66-90
Underweight	<15%	<17	<18	<19
Healthy	[15, 21]	[17, 23]	[18, 24]	[19, 25]
Overweight	[21, 26]	[23, 28]	[24, 29]	[25, 30]
Obese	>=26	>=28	>=29	>=30
Female				
Underweight	<18%	<20	<21	<22
Healthy	[18, 28]	[20, 30]	[21, 31]	[22, 32]
Overweight	[28, 35]	[30, 35]	[31, 36]	[32, 37]
Obese	>=35	>=35	>=36	>=37

Use

The device has a very sensitive sensor. If you want to adjust the settings, simply touch the buttons gently with your finger.

If the device does not respond to touch, clean the glass surface with a soft, dry cloth.

NOTE:

The scales are factory set with the units "kg" and "cm".

To change the settings to "lb" or "st" as well as "ft.in", proceed as follows:

- Switch the scales on by applying gentle pressure on the glass surface with your foot or toe.
- Press the UP button , until the required unit of weight is displayed. The unit for the height also automatically changes.

Weighing Function

- Place the scale on a hard, flat surface (avoid carpets). An uneven floor contributes to inaccurate measurements.
- Stand on the scales with both feet. The device switches on automatically.
- Distribute your weight equally and wait for the scale to determine your weight.
- The digits on the display will flash before your exact weight will be shown.
- In order to extend the battery life, the device will automatically turn off 10 seconds after you have stepped off the scale.

NOTE:

- The device can only function correctly if you stand on the device with bare, dry feet!
- If necessary wipe the contact surfaces to improve results.
- Due to the significant fluctuation of body fat, it is difficult to attain a standard value. On average the optimal testing time is between 6 p.m. and 8 p.m., since the body fat remains very consistent during this period of time.
- Press the SET button. The preset stored position data appears on the display (USER flashes). Use the UP / DOWN button to select the stored position at which your data is saved.
- Wait for the USER display and 0.0 kg (0.0 lb/st) to stop flashing.
- Stand on the contact surfaces (1) of the scales with both feet.
- Firstly your weight will flash on the display.
- After a few seconds, the following results will be shown:

1. The body fat content in % (☐ symbol)
The body water content in % (☐ symbol)
The results analysis in the form of a bar
 2. The muscle mass content in % (☐ symbol)
The bone mass content in % (☐ symbol)
The calorie requirement per day
 3. The preset gender
The preset type
The preset height
The preset age (AGE)
The weight gain (INCREASE) / the weight loss (DECREASE) in kg (lb/st) in comparison to the previous measurement
- Afterwards the determined data will be shown repeatedly. Then the scale will automatically turn off.

Display error Messages

NOTE:

If an error message appears on the display, wait until the device is turned off. Turn it on again to restart the operation. If the error display has not disappeared, take the batteries out of the battery compartment for 10 seconds.

Display Symbol	Meaning
8888	Self-test (when switching on) to check that the numeric display functions properly
----	The scale stands faultily Battery is low
Lo	Battery is low
Err	The percentage of body fat content lies above or below the limits of the measuring scale
EEEE	The weight lies above the allowed maximum weight

Cleaning

CAUTION:

- Do not immerse the balance in water; the electronic components can be damaged.
- Do not use a wire brush or any abrasive items.
- Do not use any acidic or abrasive detergents.
- Clean the device with a slightly moist cloth. Make sure that no humidity or waters seeps into the device!
- Use a commercially sold detergent or soap for cleaning.
- Dry the device well after cleaning with a dry and soft cloth.

Disposal of batteries

Your appliance contains batteries. If these batteries are spent, follow the instructions below:

CAUTION:

Batteries/rechargeable batteries must not be disposed of with household waste!

- The batteries must be removed from the appliance.
 - Take spent batteries to the appropriate collection point or to a dealer.
 - Your town or local authority can provide information about public collection points.
- This symbol can be found on batteries/rechargeable batteries which contain hazardous substances:

Pb = contains lead
Cd = contains cadmium
Hg = contains mercury
Li = contains lithium

Correct Disposal of this product

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it.

Keep this manual for future reference and pass it on with the equipment to any future users.

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Model: 756
 Power supply: 2 Batteries 1.5 V Type "AAA" / "R03"
 Net weight of device 1.3 kg
 Memory capacity: 10 PO-P9
 Mass alignment: from 100 to 220 cm
 Age: from 18 to 100 years
 Minimum weight: from 3.0 kg
 Maximum weight: 180 kg
 Weight measurement unit: 100g increments
 Body fat measurement unit: 0.1 %
 Body fat measuring range: 4.0 % - 50.0 %
 Body water measuring range: 35.0 % - 85.0 %
 Muscle mass measuring range: 10.0 % - 50.0 %
 Subject to technical changes without prior notice!

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