<b>beurer</b> <small>baby care</small> <b>BY 11</b> <span><span>CE</span></span> <span>0483</span> <span><span><span></span></span></span>	<span><span><span></span></span></span> <span>BEURER GmbH • Söflinger Str. 218 • 89077 Ulm (Germany)</span> <div>www.beurer.com</div> www.beurer.gesundheitsratgeber.com • www.beurer.healthmanager.com
DEUTSCH	

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**Important safety information – retain for later use.**

In order to make optimum use of all the merits of the clinical thermometer, you should carefully read the operating instructions before using it, retain the instructions for later use and also keep them accessible for other users.

- The thermometer is intended to measure the temperature of the human body.
- The BY 11 thermometer is only designed for the measuring area on the human body stated in the instructions for use.
- Children should not use the thermometer without supervision.
- Before each use, check the thermometer for signs of damage or wear. Damaged or worn thermometers should not be used.
- Always wait for the acoustic signal which indicates the minimum time required for measurement. Discuss the temperature with your family physician.
- The thermometer contains sensitive electronic parts. Therefore, protect it from blows, bending, high temperature or direct exposure to sunlight.
- This thermometer has a flexible measuring tip to provide greater comfort and safety during measurement, particularly in infants and persons who are asleep or have a reduced level of consciousness. The tip must not be bent more than 90°!
- Using the thermometer in strong electromagnetic fields such as next to a mobile phone, can cause a malfunction.
- When the thermometer is switched on, it carries out a self-test. Therefore the accuracy of measurement need not be tested at regular intervals.

<b>beurer</b> <small>baby care</small> <b>BY 11</b> <span><span>CE</span></span> <span>0483</span> <span><span><span></span></span></span>	<span><span><span></span></span></span> <span>BEURER GmbH • Söflinger Str. 218 • 89077 Ulm (Germany)</span> <div>www.beurer.com</div> www.beurer.gesundheitsratgeber.com • www.beurer.healthmanager.com
FRANÇAIS	

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<b>Electromagnetic Compatibility Information</b>		
<b>Table 1</b>		
<i><b>Guidance and manufacturer's declaration – electromagnetic emissions</b></i>		
The FT15/1 is intended for use in the electromagnetic environment specified below. The customer or the user of the FT15/1 should assure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The FT15/1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The FT15/1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

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**Measurement of temperature in the armpit (axillary)**

This method of measurement is relatively inaccurate. Therefore, it cannot be recommended from a medical point of view.

<b>Mode of measurement</b>	<b>Recommended time required for measurement</b>	<b>Mean divergence to the rectal method</b>
in the anus	until the acoustic signal	−
in the oral cavity	until the acoustic signal	−0.4 °C to -1.5 °C (-0.7 °F to -2.7 °F)
in the armpit	minimum 10 minutes	−0.7 °C to -2.0 °C (-1.26 °F to -3.6 °F)

### Use

Press the On/Off switch to turn the unit on. A short acoustic signal will announce that it is switched on. The thermometer will carry out a self-test for about 2 seconds. During the self-test all segments of the display will be visible. Then "Lo °C (Lo °F)" or the value of the last measurement will appear in the display. Then the reference value of 37 °C (98.6 °F) will appear. Then the flashing symbol ""C °C ("F)" appears and the symbol "Lo" is shown in the display screen. The thermometer is now ready for use. During the measurement, the current temperature is displayed continuously, and the ""C °C ("F)" symbol flashes on and off. The measurement is completed when a stable temperature has been reached. An acoustic signal (10 long beeps) sounds, the ""C °C ("F)" symbol stops flashing and the measured temperature is displayed. So long as the temperature is less than 32 °C (89.6 °F), the ""C °C ("F)" symbol continues to flash and "Lo °C (Lo °F)" is shown on the display. Higher values are displayed directly, e.g. 32.1 °C (89.8 °F). If temperatures over 37.8 °C (100 °F) are measured, a fever alarm (30 short beeps in 10 seconds) sounds. If temperatures over 43.0 °C (109.4 °F) are measured, the symbol for "Hi °C (Hi °F)" appears. The measurement time is very quick compared to conventional thermometers and is only 10 seconds for the rectal measurement. More accurate measurements are possible by extending the measurement time beyond

Press the On/Off switch to turn the unit on. A short acoustic signal will announce that it is switched on. The thermometer will carry out a self-test for about 2 seconds. During the self-test all segments of the display will be visible. Then "Lo °C (Lo °F)" or the value of the last measurement will appear in the display. Then the reference value of 37 °C (98.6 °F) will appear. Then the flashing symbol ""C °C ("F)" appears and the symbol "Lo" is shown in the display screen. The thermometer is now ready for use. During the measurement, the current temperature is displayed continuously, and the ""C °C ("F)" symbol flashes on and off. The measurement is completed when a stable temperature has been reached. An acoustic signal (10 long beeps) sounds, the ""C °C ("F)" symbol stops flashing and the measured temperature is displayed. So long as the temperature is less than 32 °C (89.6 °F), the ""C °C ("F)" symbol continues to flash and "Lo °C (Lo °F)" is shown on the display. Higher values are displayed directly, e.g. 32.1 °C (89.8 °F). If temperatures over 37.8 °C (100 °F) are measured, a fever alarm (30 short beeps in 10 seconds) sounds. If temperatures over 43.0 °C (109.4 °F) are measured, the symbol for "Hi °C (Hi °F)" appears. The measurement time is very quick compared to conventional thermometers and is only 10 seconds for the rectal measurement. More accurate measurements are possible by extending the measurement time beyond

**Measurement in the mouth cavity (oral)**

For this purpose, guide the tip of the thermometer carefully into one of the heat pouches beneath the tongue, to the left or the right of the root of the tongue.

<b>beurer</b> <small>baby care</small> <b>BY 11</b> <span><span>CE</span></span> <span>0483</span> <span><span><span></span></span></span>	<span><span><span></span></span></span> <span>BEURER GmbH • Söflinger Str. 218 • 89077 Ulm (Germany)</span> <div>www.beurer.com</div> www.beurer.gesundheitsratgeber.com • www.beurer.healthmanager.com
FRANÇAIS	

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